

Rotary evaporator RV 8

Scope of supply: Rotary evaporator RV 8, heating bath HB digital, vertical glassware RV 10.1, vapor tube, woulff bottle

Speed:	5 ... 300 1/min
Temperature:	RT ... 180 °C
Accuracy temperature:	±1 °C
Heating capacity:	1400 W
Cooling surface:	1600 cm ²
Max. volume bath:	4 l
Lifting path:	120 mm
Lift:	Manual
Condenser:	Vertical
Glassware:	Without coating
Dimensions (WxDxH):	510 x 345 x 490 mm
Weight:	25.76 kg
Power supply:	100 ... 240 V, 50/60 Hz
IP code:	IP 20

Type	Min. speed 1/min	Max. speed 1/min	Max. temperature °C	Width mm	Depth mm	Height mm	Type of auxiliary energy	PK	Cat. N
RV 8 V	5	300	180	510	345	490	Mains connection	1	9.813

Zubehör:

- Accessories for rotary evaporators RV 8/RV 10 digital
- Vacuum pump VACSTAR digital
- Vacuum pump MVP 10 basic, for rotary evaporators
- Foam brake for Rotary evaporator RV 10
- Dry ice condenser for rotary evaporator RV 10/RV 8/RV 3
- Evaporation flasks for rotary evaporator RV 10, RV 8, RV 3
- Evaporation flasks, coated for rotary evaporator RV 10/RV 8/RV 3
- Evaporation cylinders for rotary evaporator RV 10, RV 8, RV 3
- Vertical condenser with manifold and cut-off valve, for rotary evaporator RV 10/RV 8/RV 3
- Sets of glassware for Rotary evaporator RV 10
- Powder flasks for rotary evaporator RV 10, RV 8, RV 3
- Distilling spiders with distilling sleeves, for rotary evaporator RV 10, RV 8, RV 3
- Distilling spiders with flasks, for Rotary evaporator RV 10, RV 8, RV 3

Alternativen:

- Rotary Evaporator RV 8 FLEX
- Rotary evaporator RV 8 V-C